



U.S. Department of Energy
Energy Efficiency
and Renewable Energy

Bringing you a prosperous future where energy
is clean, abundant, reliable, and affordable

Surface Analysis

Sally Asher, Joel Pankow, Craig Perkins,
Robert Reedy, Glenn Teeter, and Matt Young

DOE Solar Energy Technologies Program
Peer Review

Denver, Colorado

April 17-19, 2007




Advance the understanding of photovoltaic materials and devices by employing surface science techniques to investigate surface and interfacial properties.

- **Work with Solar America Initiative (SAI) TPP awardees and other PV researchers to identify and solve research and manufacturing problems.**
- **Collaborate with DOE SETP researchers to solve key problems faced in the development and reliability of PV materials.**
- **Develop new tools and improved methodologies for *in-situ* study of PV surfaces and interfaces.**
- **Disseminate research results through publications, presentations, and regular interactions with other researchers in photovoltaic and surface-science disciplines.**



Major Analytical Methods

- Auger Electron Spectroscopy (AES)
- X-ray and Ultraviolet Photoelectron Spectroscopy (XPS/UPS)
- Surface Analysis Cluster Tool
 - UHV Transfer System
 - AES/XPS/UPS
 - UHV Deposition Tool
 - Thermal Desorption Mass Spectrometry (TDMS)
 - LEED, RHEED, SE capable
 - Controlled Ambient Glove Box
- Secondary Ion Mass Spectrometry (SIMS)
 - Dynamic (double-focusing magnetic sector)
 - Static (time-of-flight)



Phi 5600 XPS/UPS System



October 1, 2005 to January 1, 2007

| Project Task(s) | Total Value |
|--------------------|--------------------|
| PVA6.3301 | \$1,380,000 |
| PVB6.5133 | \$100,000 |
| PVA7.3301 | \$478,663* |
| Grand Total | \$1,959,663 |

FY06 - 5.88 FTE's

FY07 - 5.65 FTE's

*Amount in FY07 Plan under continuing resolution pro-rated for 1/3 year



Support Research and Development

- Worked with at least 17 different PV manufacturers
- Considerably larger number of NREL and academic researchers
- >800 samples, >3000 separate measurements
- All materials technologies

Collaborative Research and Development

- Study reaction kinetics and electronic structure of nano-scale Cu_xTe films for CdTe back contact applications
- Study relationship between moisture barrier deposition fundamentals and barrier performance
- Co-inventor of one issued US patent



Tool Development for PDIL

- **Large-area Auger tool development**
- **Liquid-Phase Quartz Crystal Microbalance (LP-QCM)**
- **Automated chemical reactor for wet-processing**
- **Theoretical calculations for flux uniformity and combinatorial deposition**
- **Large-Area XPS tool development**

Publications and Presentations


- **8 - Oral presentations - external (2 invited)**
- **2 - Oral presentations - internal**
- **12 - Publications**



Task: Study the reaction kinetics of nanoscale Cu_xTe films

Cu-doped CdTe: Temperature-Reversible Cu_xTe Surface Precipitation

Temperature-Programmed XPS Measurement



$dT/dt = 2 \text{ K/min.}$

Hysteresis
 \rightarrow nonequilibrium


XPS, AES show
 Cu_xTe at surface

Precipitation—driven
by Cu-CdTe solubility



Task: Study the electronic structure of nanoscale Cu_xTe films

Ultraviolet Photoelectron Spectroscopy and Band Alignments



| CdTe(111)-B substrate | SEC (eV) | VBM (eV) | I_s (eV) |
|-----------------------|------------|------------|------------|
| w/o Cu, T=325 K | 15.98 | 0.45 | 5.69 |
| w/ Cu, T > 635 K | 16.06 | 0.37 | 5.53 |
| w/ Cu, T=325 K | 16.38 | 0.00 | 4.84 |
| Est. uncertainty | ± 0.03 | ± 0.03 | ± 0.06 |

Heavy Cu doping does not dramatically alter valence-band or I_s .


Precipitated Cu_xTe lowers I_s by about 0.9 eV.




Task: Study the reaction kinetics of nanoscale Cu_xTe films

Ultraviolet Photoelectron Spectroscopy and Band Alignments

without Cu_xTe



with Cu_xTe





Task: Study the electronic structure of nanoscale Cu_xTe films

Ultraviolet Photoelectron Spectroscopy and Band Alignments

$\Phi_B = 0.7 \text{ eV} \rightarrow$ still too large to provide ohmic contact.


Schottky-Mott theory for MS interfaces:

$$\Phi_B = \phi_M - I_S$$

interface dipole term can decrease barrier:


$$\Phi_B = \phi_M - I_S + eV_{int}$$

Interface effects could change Φ_B









Task: Study relationship between moisture barrier deposition fundamentals and barrier performance characteristics



SiO_xN_y Moisture Barriers on PET¹

- PECVD over 12" x 12" Deposition Area
- Substrate Temperature
 - Heater setting $\neq T_{\text{Dep}}$
 - No *in-situ* annealing of PET
- Film (Non)-uniformity measured by XPS Depth Profiling
 - Thickness variation
 - Compositional variation
- *Use the film non-uniformity to perform combinatorial-type experiments for adhesion and WVTR properties.*



Still unknown what barrier properties are necessary to protect a solar cell

1. J.W. Pankow and S.H. Glick, 2006 *IEEE 4th World Conference on Photovoltaic Energy Conversion*, p. 2250.



Task: Study relationship between moisture barrier deposition fundamentals and barrier performance characteristics

SiO_xN_y Moisture Barriers on PET

- Develop Small Area WVTR MOCON Cell
 - Improved design with better thermal control and monitoring



Small Area MOCON Cell

Status as of 01/2007

- Barriers exhibit good dry adhesion (cohesive failure of test coupon rather than delamination of barrier)
 - Excellent barrier properties to ~60°C in damp heat (no detectable H₂O permeation)
 - ~>60°C lose adhesion => barrier fails
 - Performance appears independent of thickness and composition over the range studied
 - *Weak interface between barrier and PET suspected to be responsible for failure*
 - plasma damage during deposition
 - Artifact of measurement protocol
- Films are allowed to dry prior to each temperature change resulting in dimensional changes in the PET

| Sample/ Location | Barrier Composition | Thickness (Å) | T (°C) | WVTR (g/m ² day) | Peel Strength (N/mm) ¹ |
|------------------|--------------------------------------|---------------|--------|-----------------------------|-----------------------------------|
| ST504 | N/A | N/A | 40 | 3.6 | NA |
| | | | 60 | 12.5 | |
| | | | 85 | 61.0 | |
| Inlet | SiO _{1.8} N _{0.15} | 5540 | 40 | <0.05 | 0.5 |
| | | | 50 | <0.05 | |
| | | | 60 | 3.67 | |
| Center | (2) | (2) | 50 | < 0.05 | (2) |
| | | | 60 | < 0.05 | |
| | | | 70 | 8.5 | |
| Outlet | SiO _{1.6} N _{0.32} | 2755 | 50 | < 0.05 | 0.3 |
| | | | 60 | 3.7 | |

1. After 250 h damp heat
 2. To be determined
 N/A=not applicable, uncoated




Task: Design and specify a cart-based Auger electron spectroscopy tool for use in the Science and Technology Facility

Process Development and Integration Laboratory - PDIL

- Flexible platforms with integrated deposition and characterization tools
- Standard transfer interface with 6" x 6" sample size
- Rapid reconfiguration of platforms for new experiments

Large-area Auger Tool - first dedicated characterization tool

- **Whole platen analysis**
- **Highly Automated Data Acquisition**
 - Rapid depth profiling, composition
 - Stage mapping
- **Versatile chamber**
 - Ports for future experiments
- **Status as of 01/2007**
 - Plans approved
 - Under construction
 - Expected delivery 06/2007







Task: Study wet-chemical processing of PV material surfaces leading to better control of deposition processes

Liquid-Phase Quartz Crystal Microbalance (LP-QCM)

Motivation-study CBD CdS surface chemistry

- CdS thickness has strong effect on device performance (CdTe and CIGS). Thickness reproducibility critical.¹
- Post-growth thickness measurement is tedious
 - Upper right graph-Auger depth profiles from three different CdS film growths
 - Bottom right graph-three growth rate calibration curves generated from 22-separate depth profiles (represents 22 different films grown!!)
- Accuracy is poor for thinnest films
- No method for *in-situ* measurement of film growth

In-situ measurement required if PDIL tools are to be automated with endpoint detection for film growth.







Task: Study wet-chemical processing of PV material surfaces leading to better control of deposition processes
Liquid-Phase Quartz Crystal Microbalance (LP-QCM)
Simple Equipment with complex interactions in CBD

• Considerations for use

- Deposition on both sides of the crystal
 - Temperature dependent properties of system
 - Changing hydraulic pressure
 - Sensitivity to stirring
 - Bubble formation
 - Poor nucleation on sensor
 - Different reaction rates between crystal and substrate
 - Surface roughness
 - Non-rigid film
 - Film-induced stress on sensor
- ### • High sensitivity to changes in growth rate
- Ability to discern both thickness and rate during deposition
 - **R**ate is not constant throughout deposition
 - **H**owever, locally linear portions of thickness vs. time useable for endpoint detection






Task: Study wet-chemical processing of PV material surfaces leading to better control of deposition processes

Automated Reactor with LP-QCM Control of CBD Film Thickness


- **Designed a custom reactor for use in Surface Analysis cluster tool glove box**
- **Developing LabVIEW-based control software**
- **Custom hardware and software interface controls:**
 - QCM sensor
 - two chemical metering pump
 - Thermocouple
 - pH meter
 - six-solenoid mixing valve
 - spare analog inputs/outputs for expansion

Status as of 1/2007

- LP-QCM is operational
- Software controls most hardware items
- Reactor design is finalized and awaiting construction



Portion of LabVIEW .vi



Front panel for LabVIEW QCM and CBD control




Task: Perform theoretical calculations for flux uniformity and combinatorial deposition

- **Work was performed in anticipation of new large area deposition tools¹⁻²**
- **Principles were used in specification of new CIGS deposition tool**
- **Applicable to other deposition tools slated for development**

Status 1/2007

- **CIGS tool drawings have been approved**
- **Location of sources will allow combinatorial deposition per this work**



G. Teeter, *J. Vac. Sci. Technol. A*, 24(4), (2006) 1112.

G. Teeter, *J. Vac. Sci. Technol. A*, 24(4), (2006) 1119.



- **Move Surface Analysis Cluster Tool to S&TF (FY06)**
 - Completed. Tool operational September 2006
- **Large-area Auger tool for PDIL (FY06 & FY07)**
 - Tool order met FY06 milestone, Sept. 2006
 - Delivery expected June 2007, on track for FY07 milestone*
- **Large-area XPS tool for PDIL (FY07)**
 - Design specification process on track to meet FY07 milestone, Sept. 2007
- **Study wet-chemical processing of PV material surfaces leading to better control of deposition processes (FY06 & FY07)**
 - LP-QCM test and evaluation met FY06 milestone, Sept 2006
 - Progress on components of automated reactor (hardware and software components) on track to meet FY07 milestone, Aug. 2007
- **Study the reaction kinetics of nanoscale Cu_xTe films (FY06 & FY07)**
 - Completed study of Cu and Te reaction kinetics, publication 2006
 - Developing model to describe surface segregation of Cu
- **Study relationship between moisture barrier deposition fundamentals and barrier performance characteristics (FY06 & FY07)**
 - Deposition system characterized, June 2006
 - Small area MOCON cell, relationship between deposition location and barrier and adhesion properties, September 2006-ongoing
 - Deposition on Al-mirrors for barrier testing, April 2007
 - Report/project assessment due Sept. 2007
- **Investigate materials chemistry of CIGS on flexible substrates leading to an improved understanding of long-term reliability (FY07)**
 - NREL generated test samples evaluated, correlated with commercial films Nov 2006
 - Additional samples in progress
 - Report due Aug 2007



- **Continue support of, and research with SAI TPP awardees and other members of the PV community**
 - Work with M&C and NCPV management to develop prioritization and scheduling guidelines
- **Design and specify new tools for PDIL**
 - Large-area XPS tool (tool order to be placed by Sept 2007, delivery 3rd quarter FY2008)
 - Characterization cluster tool components (if funded), including wet processing station and plasma etching/thermal processing station (tools specified by Aug 2007, delivery FY2008)
- **Initiate collaborative surface science research using new large-area Auger tool**
- **Implement concepts developed from LP-QCM and automated wet-processor research into PDIL tools, specifically, semi-automated CBD CdS reactor under development by CIS task**
- **Complete quantitative model of Cu surface segregation in CdTe. Apply model to back-contact design.**
- **Determine cause for damp heat failure of barrier film interfaces. Improve coatings based on new knowledge. Recent results are promising.**
- **Investigate the effect of excess Na on the surface chemistry of CIGS, and the effect on CdS growth.**
- **Dependent on funding, upgrade surface analysis cluster tool and SIMS equipment**